

# Order form Custom Made Liner

**Fillauer**<sup>®</sup>

Kung Hans väg 2  
192 68 Sollentuna  
Sweden

Phone: +46 8 505 332 00

E-mail: fillauer@fillauer.com

Fillauer's PO#: \_\_\_\_\_

Patient name: \_\_\_\_\_

Ordered by: \_\_\_\_\_

Phone: \_\_\_\_\_



Customer: \_\_\_\_\_

## Ordering information

First order:

Repeat order:

Quantity	Type of liner	Product number First order	Product number Repeat order
	Cushion liner, black	<b>801 100 301</b>	<b>801 100 302</b>
	Cushion liner, beige	<b>801 100 401</b>	<b>801 100 402</b>
	Locking liner, black	<b>802 100 301</b>	<b>802 100 302</b>
	Locking liner, beige	<b>802 100 401</b>	<b>802 100 402</b>

Option:

	3 mm Uniform (picture 1)
	6 mm Tapered (picture 2)
	3+6 mm Contoured (picture 3)
	9 mm Contoured (only available in cushion, thick bottom, picture 4)
	Thick bottom
	Thin bottom
	Full matrix, recommended for AK prosthetic users

## IMPORTANT Information

- The plaster model needs to be a positive model.
- Make sure that the model has the same length as the final liner.
- The model should be reduced approx. 4 - 6% on circumferential measurements.
- Please state if and where you need extra cushion or other custom solutions to your liner.
- Please mark on the cast if it is a left or a right leg and the anterior/ posterior side

## Additional information:



Pic.1



Pic.2



Pic.3



Pic.4

**Fillauer**<sup>®</sup>